

## PRODUCT DATASHEET

# PILLE – Portable TL Space Dosimetry System

*4<sup>th</sup> generation*

PROJECT TITLE	-
PROJECT REF.	PRO2020
PREPARED BY	REMRED Space Technologies Ltd.
SUPERVISOR	Arisz Kecskes, kecskes.arisz@remred.hu
ISSUE	01_00
DATE OF ISSUE	01/04/2020
STATUS	Released
TYPE	Product Datasheet, Non-Confidential
REFERENCE	PRO2020-RR-PQA-DS-001_01_00
CUSTOMER(S)	-
CONTRACT REF.	-
CUSTOMER ID.	-

@REMRED 2020

The copyright in this document is vested in REMRED Ltd. This document may only be reproduced in whole or in part, or stored in a retrieval system, or transmitted in any form, or by means of electronic, mechanical, photocopying or otherwise, either with the prior permission of REMRED Ltd.

QUA2020-RR-PQA-DS-001\_01\_00

## APPROVAL

Issue	Date	Signatures		
		Author(s)	Reviews by	Approved by
01_00	01/04/2020	Andras Gerecs	Anna Baranyai	Arisz Kecskes

## DISTRIBUTION LIST

Company	Name	No. of copies
N/A	N/A	N/A
REMRED archives		1

## CHANGE LOG

Reference	Date	Issue	Revision
PRO2020-RR-PQA-DS-0001	01/04/2020	01_00	1

## CHANGE RECORD

Issue: 01_00	Revision: 1		
Reason for change	Date	Page	Paragraph(s)
N/A (initial release).	01/04/2020	All	All

## TABLE OF CONTENTS

<b>1</b>	<b>Purpose and Scope</b> .....	<b>4</b>
<b>2</b>	<b>Application and Key Features</b> .....	<b>5</b>
	2.1 APPLICATION .....	5
	2.2 KEY FEATURES.....	5
<b>3</b>	<b>Specification</b> .....	<b>6</b>
	3.1 GENERAL SPECIFICATION .....	6
	3.2 ENVIRONMENTAL SPECIFICATION .....	6
	3.3 INTERFACES .....	7
	3.4 MEASUREMENT CAPABILITIES.....	7
<b>4</b>	<b>Flight Heritage</b> .....	<b>8</b>
<b>5</b>	<b>List of Abbreviations</b> .....	<b>9</b>
<b>6</b>	<b>List of Figures</b> .....	<b>10</b>
<b>7</b>	<b>List of Tables</b> .....	<b>11</b>
<b>8</b>	<b>References</b> .....	<b>12</b>
	8.1 APPLICABLE AND NORMATIVE DOCUMENTS.....	12
	8.2 REFERENCE DOCUMENTS.....	12

# 1 Purpose and Scope

The present document provides detailed technical information about the PILLE Thermoluminescent Space Dosimetry System.

The definitions and glossary of terms from ECSS-S-ST-00-01C [AD 1] apply to this document.



Figure 1 – PILLE System (4<sup>th</sup> generation)

©REMRED 2020

The copyright in this document is vested in REMRED Ltd. This document may only be reproduced in whole or in part, or stored in a retrieval system, or transmitted in any form, or by means of electronic, mechanical, photocopying or otherwise, either with the prior permission of REMRED Ltd.

QUA2020-RR-PQA-DS-001\_01\_00

## 2 Application and Key Features

### 2.1 APPLICATION

- ✓ **Space dosimetry monitoring during manned space missions**
  - A unique, proven space dosimetry instrument for manned space missions
  - To measure the absorbed dose at different locations in the space vehicle/station (space dosimetry mapping)
  - To study radiation shielding effects of the surrounding environment
  - To provide early warnings for astronauts about dose levels
  - Utilized for routine and extravehicular activity (EVA) individual dosimetry of astronauts

### 2.2 KEY FEATURES

- ✓ **Space dosimetry instrument for manned space missions to determine the absorbed dose from space radiation**
- ✓ **High reliability in space environment**
  - Proven >10years operation in space (on-board ISS at LEO)
- ✓ **Thermoluminescent detector system (TLD) including detector(s) and reader**
- ✓ **Unlimited number of passive detectors with memory chip inside containing the identification code and the individual calibration parameters of the dosimeter**
- ✓ **Easy-to-use, compact, lightweight microprocessor controlled on-board Reader Unit for providing the preliminary evaluation of the dose absorbed by the dosimeters**
- ✓ **The Reader Unit provides**
  - User interface for astronauts (via numerical display and keyboard)
  - Automatic identification of the detectors
  - The measured dose and a series of parameters are displayed and stored on a removable memory card
  - A dosimeter inserted in the reader permanently is dedicated for automatic cyclic measurements
- ✓ **RS-232 interface to connect to a PC for using a dedicated PC software (“PILLE Controller”) for downloading, listing and displaying the results of the measurements and setting the parameters of the Reader and the dosimeters**
- ✓ **CAN interface to connect to the data acquisition system of the space station**
- ✓ **Operated on-board space stations Salyut-6,-7; Mir and the ISS**

# 3 Specification

## 3.1 GENERAL SPECIFICATION

1. Table – General specification

Parameters	Reader	Detector
Power	0.1/1/7W (standby/ready/readout)	-
Mass	~ 1.4 kg	70g (with carrying case)
Dimensions (H, W, D)	70mm, 190 mm, 120 mm	ø 20 mm * 60 mm
Operational temperature range	-20°C...+40°C	-40°C...+50°C
Non-operatioinal temperature range	-40°C...+85°C	-40°C...+85°C
Operational pressure range	1.2 10 <sup>5</sup> Pa...7 10 <sup>-4</sup> Pa	2 10 <sup>5</sup> Pa... <sup>-4</sup> Pa
Data rate	512 byte / readout	-

## 3.2 ENVIRONMENTAL SPECIFICATION

2. Table – Environmental specification

Parameters	Values	
Thermal	Temperature environment	-40°C...+85°C
	Pressure environment	Up to 1.2 10 <sup>5</sup> Pa
	Humidity environment	25...65%relH
Thermal-vacuum	Temperature environment	-40°C...+85°C
	Vacuum environment	<10 <sup>-3</sup> Pa
	Max. depressurisation rate	5.0 kPa/s
Vibration	Sine vibration environment for 3-axis	20...100 Hz, 16.0 g
	Random vibration environment for 3-axis	5...2000Hz, 17.0g <sup>RMS</sup>
	Shock pulse for 3 axis	100 g, 0.25 ms
EMC	EMC environment	Tailored ECSS-E-ST-20-07C Rev.1 [AD 2]
Radiation	Used components	COTS
	Proven lifetime	>10 years proven lifetime in LEO

## 3.3 INTERFACES

3. Table – Interfaces

Parameters	Values
Input power bus*	17.0 V...34V (non-redundant)
TM/TC interface*	CAN Bus / RS-232 (redundant/parallel)
*For Reader only	

## 3.4 MEASUREMENT CAPABILITIES

4. Table – Measurement capabilities

Dosimeters	
Type	bulb
Material	CaSO <sub>4</sub> :Dy
Reader	
Measuring range (s<10%)	3μGy...10 Gy
TLD Efficiency (e=1±%)	LET <sub>∞H20</sub> < 10 keV/μm
Read-out precision	3 digits + exp.
Accuracy (above 10 μGy)	δ < 5%
Measuring modes	manual / automatic read-out
Display	8-digit alphanumeric LED
Displayed information	dose in μGy (air kerma); date and time of measurement; identification codes; mean dose rate; error codes
Storage of information	PCMCIA memory card (> 4000 data sets)

## 4 Flight Heritage

5. Table – Flight heritage

Mission name	Hosting platform	Orbit details	Duration	Remarks
PILLE	Salyut-6 SS	LEO 200-300 km	1,5 years	16 dosimeters and one Reader
PILLE	Salyut-7 SS	LEO 300-400 km	3 years	16 dosimeters and one Reader
PILLE'S	STS 41G / Challenger-6	LEO 300-400 km	8 years	8 dosimeters and one Reader
PILLE	Mir SS	LEO 300-400 km	4 years	16 dosimeters and one Reader
PILLE'95	Mir SS	LEO 300-400 km	5 months	6 dosimeters and one Reader
PILLE'96	Mir SS (NASAMir4)	LEO 300-400 km	5 months	8 dosimeters and one Reader
PILLE ISS	ISS – Destiny (U.S. Lab)	LEO 330-420 km	3 month / 10 years	50 dosimeters and one Reader
PILLE-MKS	ISS – Zvezda (Rus.Lab)	LEO 330-420 km	15 years	12-17 dosimeters and one Reader



## 5 List of Abbreviations

<b>AD</b>	Applicable Documents
<b>COTS</b>	Commercial Off-The-Shelf
<b>ECSS</b>	European Cooperation for Space Standardization
<b>EVA</b>	ExtraVehicular Activity
<b>ISS</b>	International Space Station
<b>LEO</b>	Low Earth Orbit
<b>LET</b>	Linear Energy Transfer
<b>PC</b>	Computer
<b>RD</b>	Reference Documents
<b>SS</b>	Space Station
<b>LT</b>	Thermoluminescent
<b>TLD</b>	Thermoluminescent detector system

# 6 List of Figures

**Figure 1 – PILLE System (4<sup>th</sup> generation) .....4**

## 7 List of Tables

1. Table – General specification .....	6
2. Table – Environmental specification .....	6
3. Table – Interfaces .....	7
4. Table – Measurement capabilities .....	7
5. Table – Flight heritage .....	8
6. Table – Applicable and Normative Documents .....	12
7. Table – Reference Documents .....	12

# 8 References

## 8.1 APPLICABLE AND NORMATIVE DOCUMENTS

6. Table – Applicable and Normative Documents

AD	Title	Reference	Issue
[AD 1]	ECSS system - Glossary of terms	ECSS-S-ST-00-01C	1 Oct 2012
[AD 2]	Space engineering – Electromagnetic compatibility	ECSS-E-ST-20-07C Rev.1	7 Feb 2012

## 8.2 REFERENCE DOCUMENTS

7. Table – Reference Documents

RD	Title	Reference	Issue
[RD 1]	-	-	-

**END OF DOCUMENT**